


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED



ITEM	QTY.	PART NO.	REV	CODE	NAME	DESCRIPTION
1	1	80664	01	UPSTRY	CHEST VEST SMALL	STRAP CHEST VEST SM
2	4	81148	01	UPSTRY	CHEST VEST STRAP	CHEST VEST STRAP

NOTES:

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES AND SPECS. ARE:	CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE DO NOT SCALE DRAWING	ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900	
DECIMALS .XXX ± .031	ANGLES ± .5°	APPROVALS	DATE
BREAK EDGES .005 - .015	HOLE SIZE TOLERANCE: (+.005 AND -.000)	DRAWN: MWL	03/03/03
ADD A .031 DEEP 45° CHAMFER ON ALL DRILLED HOLES	CODE: UPSTRY	MFG MNG:	MFG ENG:
FINISH: -	MATERIAL: -	ASMBLY:	SHIPPING:

NAME: STRAP CHEST VEST SM EI	REV. 000
DESCRIPTION: STRAP CHEST VEST SM EI	
SIZE DWG. NO. 81149	
SCALE = 1:3	CAD FILE NAME: 81149
SHEET 1	OF 1