


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED



NOTES:		CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE		ALTIMATE MEDICAL, INC.	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES AND SPECS. ARE:		DO NOT SCALE DRAWING		PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900	
DECIMALS .XXX ± .031	ANGLES ± .5°	APPROVALS	DATE	NAME:	
BREAK EDGES .005 - .015	HOLE SIZE TOLERANCE: (+.005 AND -.000)	DRAWN: MWL	03/03/03	STRAP CHEST VEST MED EI	
ADD A .031 DEEP 45° CHAMFER ON ALL DRILLED HOLES		MFG MNG:	MFG ENG:	DESCRIPTION:	
CODE: UPSTRY	MATERIAL: NYLON	ASMBLY:	SHIPPING:	SIZE DWG. NO. B	81150
FINISH: -	SCALE = 1:3		CAD FILE NAME: 81150	REV. 000	SHEET 1 OF 1

ITEM	QTY.	PART NO.	REV	CODE	NAME	DESCRIPTION
1	1	80665	A	UPSTRY	CHEST VEST MEDIUM	STRAP CHEST VEST MED
2	4	81148	01	UPSTRY	CHEST VEST STRAP	CHEST VEST STRAP