


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	05-ECO-001	7/12/05	JVS



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG30019	000	1	VEST CHEST Y SHAPE LARGE
2	NG30020	A	4	BUCKLE SNAP W/LOOP

- IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.
- SEE SHIPPING INSTRUCTIONS
 - RECEIVING INSPECTION
 - SEE ENVIRONMENTAL INSTRUCTIONS

APPROVED:		DATE:
DRAWN:	DATE:	MATERIAL:
JVS	3/2/2005	BLK
OPTION NUMBER:	CODE:	
V903	UPSTRY	

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:

DECIMALS	ANGLES
.X ± .1	.X ± .5°
.XX ± .03	
.XXX ± .015	
HOLE SIZE TOLERANCE: (+.005/- .000)	
BREAK EDGES .005 - .015	
CMFR DRILLED HOLES .03 X 45°	

ALTIMATE MEDICAL, INC.
 PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270
 PHONE (507) 697-6393 FAX (507) 697-6900

DESCRIPTION:
VEST CHEST Y SHAPE ASSEMBLY LARGE

SIZE DWG. NO. **NG30021** REV. **000**

SCALE = 1:4. \$CAD FILE NAME: NG30021 SHEET 1 OF 1