


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	05-ECO-001	8/15/05	GJR
A	REVISED ITEM NG30022.	DCR-00089	11/27/06	MWL
B	REVISED ITEM NG30022.	DCR-00155	3/16/07	MWL



NOTES:
 1. ATTACH ROHO CUSHION TO SEAT USING VELCRO TABS PROVIDED FROM ROHO.

ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG30013	000	1	PAD SEAT ROHO CUSHION
2	NG30022	B	1	PAD SEAT LG ROHO RECEIVER

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE: DECIMALS ANGLES .X ± .1 .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900 DESCRIPTION: PAD SEAT LG ROHO OPTION SIZE DWG. NO. NG30023 REV. B SCALE = 1:6 CAD FILE NAME: NG30023 SHEET 1 OF 1	
	DRAWN: MWL	DATE: 8/15/05			MATERIAL:
	OPTION NUMBER: V705	FINISH: -			CODE: UPSTRY