


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	06-ECO-051	11/27/06	MWL
A	REVISED NG30055.	DCR-00155	3/16/07	MWL



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG30055	A	1	PAD SEAT ROHO YOUTH
2	81035	-	1	ROHO CUSHION

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	<b>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:</b> DECIMALS      ANGLES .X ± .1      .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	<b>ALTIMATE MEDICAL, INC.</b> PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393      FAX (507) 697-6900 DESCRIPTION: <b>PAD SEAT ROHO YOUTH</b> SIZE DWG. NO. <b>NG30183</b> REV. <b>A</b> SCALE = 1:3      CAD FILE NAME: NG30183      SHEET 1 OF 1	
	DRAWN: JVS	DATE: 9/18/2006			MATERIAL:
	OPTION NUMBER: -	FINISH: -			CODE: ASMBLY