


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE	DCR-00804	6/24/11	TJF



GENERAL MATERIAL PROPERTIES

EXTERIOR COVERING PROPERTIES:

1. MATERIAL: DARTEX P091, P191 AND P291 (C, S & P) TYPES
2. DESCRIPTION: POLYURETHANE TRANSFER COATING ON WEFT KNITTED FABRIC
3. MATT GLOSS LEVLE
4. COATING PROTECTED AGAINST FUNGI AND BACTERIA
5. WASHABLE
6. DRY CLEANABLE
7. HF WELDABLE
8. FLAMABILITY: BS 7175 IGNITION SOURCES 0, 1 & 5: CONFORMS TO CALIFORNIA 117 SECTION E: COMPLIES (NON-IGNITION)

NEOPRENE PROPERTIES:

9. 6mm THICK NYLON FABRIC-FACED NEOPRENE
10. HARDNESS DUROMETER: 61-65 SHORE A
11. TENSILE STRENGTH - 85 psi
12. STRETCH LIMIT - 450%
13. DENSITY 11 lbs/cubic foot
14. FLAME RESISTANCE TO MEET CAL 117

NOTES:

1. FEED DRAWSTRING THROUGH THE SLEEVE ON THE BOTTOM OF THE TRAY COVER AND THROUGH BOTH DRAWSTRING STOPPERS.
2. DRAWSTRING STOPPERS TO KEEP THE TRAY COVER TIGHT ON THE TRAY

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECIEVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE: DECIMALS ANGLES .X ± .1 .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900		
	DRAWN: TJF	DATE: 6/22/11		MATERIAL: DARTEX	DESCRIPTION: PADDED TRAY COVER	
	EXPIRATION DATE:			FINISH: -	CODE: UPSTRY	SIZE DWG. NO. NG30295
			SCALE = 1:4	CAD FILE NAME: NG30295	SHEET 1 OF 1	