



THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
01	ADDED EXTERNAL TOOTH (STAR) WASHERS	-	4/20/05	GJR
000	NEW RELEASE.	05-ECO-001	7/12/05	JVS
-	REMOVED STAR WASHERS #600141	-	11/8/05	JVS
A	REVISED #600005 TO REV B	DCR-00699	8/24/10	JVS
B	ADDED NG40058/DECAL	DCR-00919	5/7/12	JVS



ITEM NO.	PART NO.	REV	QTY.	NAME
1	600005	B	2	CASTER 125mm BLACK
2	NG10019	000	1	BRACKET CASTER MOUNT LG
3	NG40058	000	1	DECAL "SEE INSTRUCTIONS FOR USE"

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	<b>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:</b> DECIMALS .X ± .1    ANGLES .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/- .000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	<b>ALTIMATE MEDICAL, INC.</b> PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393    FAX (507) 697-6900 DESCRIPTION: <b>SEAT REAR CASTER TUBE</b> SIZE DWG. NO. <b>NG50058</b> REV. <b>B</b> SCALE = 1:4    CAD FILE NAME: NG50058    SHEET 1 OF 1	
	DRAWN: JVS	DATE: 4/6/2005			MATERIAL:
					FINISH: -