


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	05-ECO-001	8/25/05	MWL
A	REVISED ITEM NG10097.	DCR-00161	3/20/07	MWL
B	REVISED SOCKET ANGLE ON NG10097.	DCR-00471	12/12/08	MWL
C	REVISED NG10097 TO REV C, REPLACED 91501 WITH 91406 & 91511	DCR-00570	7/29/09	JVS



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG10097	C	1	PUMP HANDLE WELDMENT RT
2	91406	000	2	BOLT 1/4-20 X 1-1/2 SHCS
3	91511	000	2	NUT 1/4-20 NYLON LOCK

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE: DECIMALS ANGLES .X ± .1 .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900 DESCRIPTION: HANDLE PUMP STANDARD RT SIZE DWG. NO. NG50059 REV. C SCALE = 1:6 CAD FILE NAME: NG50059 SHEET 1 OF 1	
	DRAWN: JVS	DATE: 4/6/2005			MATERIAL:
					CODE: ASMBLY