


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

| REVISIONS |                                |            |          |           |
|-----------|--------------------------------|------------|----------|-----------|
| REV:      | DESCRIPTION:                   | ECO #:     | DATE:    | APPROVED: |
| 000       | NEW RELEASE.                   | 05-ECO-001 | 8/23/05  | GJR       |
| A         | ADDED HOLE PLUGS 600255        | 06-ECO-029 | 5-17-06  | GJR       |
| B         | REPLACED NG50086 WITH NG50165. | 06-ECO-051 | 11/21/06 | MWL       |
| C         | REVISED NG50165 TO REVISION A. | DCR-00111  | 1/11/07  | MWL       |
| D         | REVISED NG50165 TO REVISION B. | DCR-00243  | 8/6/07   | MWL       |
| E         | REVISED NG50165 TO REVISION C. | DCR-00291  | 11/15/07 | MWL       |
| F         | REVISED NG50165 TO REVISION D. | DCR-00395  | 4/23/08  | MWL       |
| G         | REORGANIZED BILL OF MATERIALS. | DCR-00459  | 10/27/08 | MWL       |



| ITEM NO. | PART NO. | REV | QTY. | NAME                         |
|----------|----------|-----|------|------------------------------|
| 1        | NG30049  | 000 | 1    | PAD BACK REST CONTOURED HIGH |
| 2        | 600060   | 000 | 3    | BOLT M8 X 1.25 X 50 BTN      |
| 3        | NG10057  | A   | 1    | TUBE BACK REST HIGH          |
| 4        | 600255   | 000 | 2    | CAP END M8 FLUSH             |
| 5        | 600037   | 000 | 2    | BOLT M6 X 1.0 - 10 BTN SHCS  |

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.

SEE SHIPPING INSTRUCTIONS

RECEIVING INSPECTION

SEE ENVIRONMENTAL INSTRUCTIONS

APPROVED: \_\_\_\_\_ DATE: \_\_\_\_\_

DRAWN: JVS DATE: 4/26/2005 MATERIAL: \_\_\_\_\_

FINISH: \_\_\_\_\_

EXPIRATION DATE: \_\_\_\_\_ CODE: ASMBLY

UNLESS OTHERWISE SPECIFIED  
 DIMENSIONS ARE IN INCHES.  
 TOLERANCES AND SPECS. ARE:

DECIMALS ANGLES  
 .X ± .1 .X ± .5°  
 .XX ± .03  
 .XXX ± .015

HOLE SIZE TOLERANCE:  
 (+.005/-.000)

BREAK EDGES .005 - .015  
 CMFR DRILLED HOLES .03 X 45°

**ALTIMATE MEDICAL, INC.**  
 PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270  
 PHONE (507) 697-6393 FAX (507) 697-6900

DESCRIPTION:  
 BACK HIGH W/CONTOUR (V910)

SIZE DWG. NO. **NG50068** REV. **G**

SCALE = 1:6 CAD FILE NAME: NG50068 SHEET 1 OF 1