


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	06-ECO-051	9/15/06	MWL



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG10727	000	2	HANDLE ADJUSTMENT VERTICAL MEDIUM
2	NG10726	000	2	WLDMT GLIDER MED HANDLE
3	92703	00	2	BALL CAP
4	600220	000	4	BOLT BTN SHCS M6 X 1.0 - 12
5	90101	00	1	24" FOAM GRIP
6	90101	00	1	24" FOAM GRIP

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.		APPROVED:		DATE:		UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:		<b>ALTIMATE MEDICAL, INC.</b> PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900	
<input type="checkbox"/> SEE SHIPPING INSTRUCTIONS	<input type="checkbox"/> RECIEVING INSPECTION	<input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	DRAWN: JVS	DATE: 10/18/2006	MATERIAL:	FINISH:	CODE: ASMBLY	DECIMALS: .X ± .1 ANGLES: .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/- .000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	DESCRIPTION: <b>HANDLE ADJ U-SHAPE YOUTH</b>
OPTION NUMBER: -								SIZE DWG. NO. <b>NG50251</b> REV. 000 SCALE = 1:5 CAD FILE NAME: NG50251 SHEET 1 OF 1	