

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE	DCR-00475	12/23/08	JVS
A	REVISED NG50454 TO REVISION C.	DCR-00827	9/8/11	TJF
B	REVISED NG50454 TO REVISION D	DCR-00905	3/14/12	JVS



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG11086	000	1	TABLE H-FRAME EVOLV
2	80825	000	1	SNAP BUTTON
3	NG50454	D	1	TRAY ASSEMBLY BLK-MOLDED H-FRAME

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	<b>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:</b> DECIMALS      ANGLES .X ± .1          .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	
	DRAWN: JVS	DATE: 11/19/08		MATERIAL:
	EXPIRATION DATE: -	FINISH: -		CODE: ASMBLY
ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393      FAX (507) 697-6900		DESCRIPTION: <b>TABLE ASSEMBLY BLK MOLDED H-FRAME</b>		
SIZE DWG. NO. <b>NG50418</b>		REV. <b>B</b>		
SCALE = 1:5    CAD FILE NAME: NG50418		SHEET 1 OF 1		