

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE	DCR-00475	12/23/08	JVS
A	REVISED NG50457 TO REV C	DCR-00905	3/15/12	JVS



ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG11086	000	1	TABLE H-FRAME EVOLV
2	80825	000	1	SNAP BUTTON
3	NG50457	C	1	TRAY ASSEMBLY OVERSIZED H-FRAME

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.

- SEE SHIPPING INSTRUCTIONS
- RECEIVING INSPECTION
- SEE ENVIRONMENTAL INSTRUCTIONS

APPROVED: _____ DATE: _____

DRAWN: JVS DATE: 11/21/08 MATERIAL: _____
 FINISH: _____
 CODE: ASMBLY
 EXPIRATION DATE: -

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES.
 TOLERANCES AND SPECS. ARE:

DECIMALS .X ± .1 ANGLES .X ± .5°
 .XX ± .03
 .XXX ± .015
 HOLE SIZE TOLERANCE:
 (+.005/-.000)
 BREAK EDGES .005 - .015
 CMFR DRILLED HOLES .03 X 45°

ALTIMATE MEDICAL, INC.
 PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270
 PHONE (507) 697-6393 FAX (507) 697-6900

DESCRIPTION:
 TABLE ASSEMBLY OVERSIZED H-FRAME

SIZE DWG. NO. **NG50458** REV. **A**
 SCALE = 1:5 CAD FILE NAME: NG50458 SHEET 1 OF 1