

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	DCR-00593	1/27/10	MWL
A	REPLACED BOLT 600197 WITH 600022.	DCR-00658	3/26/10	MWL



ITEM NO.	PART NO.	REV	QTY.	NAME
1	T10285	000	1	BRACKET HIGH MOUTN CHEST VEST
2	600022	000	4	SHCS M6 X 1.0-20
3	600008	000	4	NUT M6 X 1.0 NYLON LOCK

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.		APPROVED:	DATE:
<input type="checkbox"/> SEE SHIPPING INSTRUCTIONS	<input type="checkbox"/> RECIEVING INSPECTION	DRAWN: MWL	DATE: 9/2/09
<input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS		EXPIRATION DATE:	-
		MATERIAL:	
		FINISH:	-
		CODE:	ASMBLY

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES.
 TOLERANCES AND SPECS. ARE:

DECIMALS	ANGLES
.X ± .1	.X ± .5°
.XX ± .03	
.XXX ± .015	

HOLE SIZE TOLERANCE:
 (+.005/-.000)

BREAK EDGES .005 - .015
 CMFR DRILLED HOLES .03 X 45°

ALTIMATE MEDICAL, INC.
 PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270
 PHONE (507) 697-6393 FAX (507) 697-6900

DESCRIPTION:
 ASMBLY HIGH MOUNT CHEST VEST BRKT

SIZE DWG. NO. **T50128** REV. **A**

SCALE = 1:4 CAD FILE NAME: T50128 SHEET 1 OF 1