


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE	DCR-00729	11/9/10	TJF
A	REVISED T50152.	DCR-00740	12/8/10	MWL
B	REVISED T50152	DCR-00807	7/7/2011	JVS
-	NOTE: T50152 TO REV C WAS TO UPDATE STICKERS	-	-	-
C	REVISED T50152 TO REV D	DCR-00925	5/29/12	TJF



**NOTE:**  
 1. MUST PULL TWO (2) WARNING LABELS WRITTEN EITHER IN ENGLISH(NG40033) OR FRENCH (NG40048) AND PLACE ONE ON EACH SIDE AS SHOWN.

ITEM NO.	PART NO.	REV	QTY.	NAME
1	T10270	000	1	FRAME MOBILE BANTAM
2	T50152	D	2	ASSEMBLY MOBILE BANTAM

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT. <input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	APPROVED:	DATE:	<b>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:</b> DECIMALS      ANGLES .X ± .1          .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	<b>ALTIMATE MEDICAL, INC.</b> PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393      FAX (507) 697-6900		
	DRAWN: MWL	DATE: 6/14/10		MATERIAL:	DESCRIPTION: MOBILE OPTION BANTAM	
	EXPIRATION DATE: -	FINISH: -		CODE: ASMBLY	SIZE DWG. NO. T50150	REV. C
				SCALE = 1:5	CAD FILE NAME: T50150	SHEET 1 OF 1