


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT  
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE	DCR-00806	7/6/2011	JVS
-	NOTE: T50178 TO REV A WAS TO UPDATE STICKERS	-	-	-
A	REVISED T50178 TO REV B	DCR-00925	5/29/12	TJF



**NOTE:**  
 1. MUST PULL TWO (2) WARNING LABELS WRITTEN EITHER IN ENGLISH(NG40033) OR FRENCH (NG40048) AND PLACE ONE ON EACH SIDE AS SHOWN.

ITEM NO.	PART NO.	REV	QTY.	NAME
1	T10270	000	1	FRAME MOBILE BANTAM
2	T50178	B	2	ASSEMBLY MOBILE BANTAM TALL

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.	APPROVED:	DATE:	UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE: DECIMALS .X ± .1 ANGLES .X ± .5° .XX ± .03 .XXX ± .015 HOLE SIZE TOLERANCE: (+.005/-.000) BREAK EDGES .005 - .015 CMFR DRILLED HOLES .03 X 45°	<b>ALTIMATE MEDICAL, INC.</b> PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900 DESCRIPTION: MOBILE OPTION BANTAM TALL SIZE DWG. NO. <b>T50177</b> REV. <b>A</b> SCALE = 1:5 CAD FILE NAME: T50177 SHEET 1 OF 1
<input type="checkbox"/> SEE SHIPPING INSTRUCTIONS <input type="checkbox"/> RECEIVING INSPECTION <input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS	DRAWN: JVS	DATE: 7/5/2011		
	EXPIRATION DATE:	-	FINISH:	
			CODE: ASMBLY	