


THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF
 ALTIMATE MEDICAL. ANY REPRODUCTION IN PART OR WHOLE WITHOUT
 THE WRITTEN PERMISSION OF ALTIMATE MEDICAL IS PROHIBITED.

REVISIONS				
REV:	DESCRIPTION:	ECO #:	DATE:	APPROVED:
000	NEW RELEASE.	DCR-01002	11/19/12	MWL



REFERENCE Y5594 OR Y5596

ITEM NO.	PART NO.	REV	QTY.	NAME
1	NG20134	000	2	KNEE COVER BACK INDEPENDENT
2	NG30290R	000	1	KNEE ROHO PAD ASSY IND RIGHT
3	NG30290L	000	1	KNEE ROHO PAD ASSY IND LEFT

IF BOX CONTAINS AN "X", SEE NOTES THAT APPLY TO EACH SUBJECT.		APPROVED: _____ DATE: _____		UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES AND SPECS. ARE:		ALTIMATE MEDICAL, INC. PO BOX 180 - 262 WEST FIRST STREET - MORTON, MN 56270 PHONE (507) 697-6393 FAX (507) 697-6900	
<input type="checkbox"/> SEE SHIPPING INSTRUCTIONS	<input type="checkbox"/> RECEIVING INSPECTION	DRAWN: MWL	DATE: 5/29/12	MATERIAL:	DECIMALS .X ± .1	ANGLES .X ± .5°	DESCRIPTION: ASSEMBLY KNEE PADS ROHO
<input type="checkbox"/> SEE ENVIRONMENTAL INSTRUCTIONS		EXPIRATION DATE: -		FINISH: -	HOLE SIZE TOLERANCE: (+.005/-.000)		SIZE DWG. NO. Y5600
				CODE: ASMBLY	BREAK EDGES .005 - .015		REV. 000
					CMFR DRILLED HOLES .03 X 45°		SCALE = 1:4 CAD FILE NAME: Y5600 SHEET 1 OF 1