

Bill TO (Ultimate Medical Supplier Only)

Purchase Order #:
Quote #:
Account #:
Ordered By:
Supplier Company:
Address:
City, State, Zip:
Phone: _____ Email: _____
Payer Source:
Name of ATP (Required):

Ship To

Name:
Address:
City, State, Zip:
Phone:
Tag For:

Bantam Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat pivot to planar back)	15"-20"
Seat to Foot Plate Range (from seat pivot)	7"-18.5"

Base-must select one Note: The DEFAULT items will be added if no selection is made.

PY5690 EasyStand[®] Bantam Medium Supine Ready Base \$2,274
 -Supine CANNOT be ordered off-unit

Supine-must select one

- PY5688 Supine (head support required) \$559
- PY5689 Non-Supine N/C

PY5691 EasyStand[®] Bantam Medium Base \$2,274
 -Supine not available with PY5691

Choose actuator handle location

- PY5698 Right Side DEFAULT N/C
- PY5699 Left Side N/C

Lift Mechanisms Accessories

- PY5632 Adjustable/Removable Actuator Handle (not available with PY5690) \$236
- PNG50044 Quad Grip Handle Extension (not available with PY5690) \$159
- PNG50070 T-Style Handle Extension (not available with PY5690) \$159

Front Frame-must select one

- PY5694 Front Frame (5" Front Wheels) **DEFAULT** \$163
- PY5695 Front Frame (3" Front Swivel Casters) \$204

Shadow Trays-must select one

- PY5604 Black Molded Tray-21"Lx24"W (depth, height and angle adjustable) \$686
- PY5606 Clear Tray-21"Lx24"W (depth, height and angle adjustable) \$948
- PY5656 Black Molded Swing-Away Tray-21"Lx24"W **DEFAULT** \$813
(depth, height and angle adjustable)
- PY5657 Clear Swing-Away Tray-21"Lx24"W (depth, height and angle adjustable) \$1,003
- PNG30295 Padded Tray Cover (not available with elbow pads) \$234

Tray Accessories

- PNG50042 Hand Grips (for clear trays only, pair) \$164
- PT50156 Elbow Pad-9.25" Cutout \$155
- PT50160 Elbow Pad-11" Cutout \$155
- PY5582 Elbow Pad-13" Cutout \$155
- PNG30000 Large Contoured Chest Pad (not available with elbow pads) \$147

Knee Pads-must select one

- PY5594-1 Slide-In Knee Pads (pair) \$435
- PY5596-1 Swing-Away Knee Pads (pair) **DEFAULT** \$701

Knee Pad Size-must select one

- PY5624 4.0" N/C
- PY5598 5.5" **DEFAULT** N/C
- PY5600-1 5.5" Roho[®] Knee Pads (pair) \$323

Foot Plates-must select one

- PY5584-1 Height-Adjustable Foot Plates (pair) \$308
- PY5586-1 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$342

Foot Plate Size-must select one

- PY5588 7.75"Lx3.25"W (pair) N/C
- PY5590 9.75"Lx4"W (pair) N/C
- PY5592 11.75"Lx5"W (pair) **DEFAULT** N/C

Foot Plate Accessories

- PT50080 Foot Straps-10"L (length over top of foot, pair) \$50
- PY5636 Foot Straps-15"L (length over top of foot, pair) \$50
- PT50082 Secure Foot Straps-10"L (length over top of foot, two pair) \$100
- PY5638 Secure Foot Straps-15"L (length over top of foot, two pair) \$100

Seats-must select one

- PY5522-1 Planar Seat **DEFAULT** \$402
- PY3042 Black Hygienic Seat Cover for PY5522-1 \$127
- PY5570 Contoured Seat-Narrow 11"-13"W \$402
- PY3043 Black Hygienic Seat Cover for PY5570 \$127
- PY5572 Contoured Seat-Wide 12"-14"W \$402
- PY3044 Black Hygienic Seat Cover for PY5572 \$127

Seat Accessories

- PY5645 Easy-Adjust Seat Depth-CANNOT be ordered off-unit \$112
(roller bearing allows for smoother seat depth adjustment)

Pelvic Support

<input type="checkbox"/>	PY5560	Hip Supports (9"-15" range, pair, not available with contoured seat)	\$376
<input type="checkbox"/>	PY3008	Velcro® Positioning Belt (fits hip circumference up to 66", velcro with D-ring)	\$65
<input type="checkbox"/>	PY3014	Positioning Belt with Airline Style Buckle (fits hip circumference up to 66")	\$159

Backs-must select one

<input type="checkbox"/>	PY5530	Planar Back-16"H	\$757
<input type="checkbox"/>	PY3041	Black Hygienic Back Cover for PY5530	\$94
<input type="checkbox"/>	PY5532	Planar Back-20"H DEFAULT	\$800
<input type="checkbox"/>	PY3040	Black Hygienic Back Cover for PY5532	\$94
<input type="checkbox"/>	PY5574	Contoured Back-16"H	\$835
<input type="checkbox"/>	PY3045	Black Hygienic Back Cover for PY5574	\$94
<input type="checkbox"/>	PY5576	Contoured Back-20"H	\$877
<input type="checkbox"/>	PY3046	Black Hygienic Back Cover for PY5576	\$94

Back Accessories

<input type="checkbox"/>	PNG50040	Push Handles (pair)	\$159
<input type="checkbox"/>	PY5643	Angle Locator	\$40

Lateral Supports

<input type="checkbox"/>	PY5686	Lateral Supports with Flat Pads-9"-14.5" Range (for planar back only, pair)	\$376
<input type="checkbox"/>	PY5562	Lateral Supports with Curved Pads-9"-14.5" Range (for planar back only, pair)	\$376
<input type="checkbox"/>	PY5658	Lateral Supports with Elbow Stop-9"-14.5" Range (for planar back only, pair)	\$384
<input type="checkbox"/>	PY5659	Lateral Supports with Elbow Stop and Arm Rest-9"-14.5" Range (for planar back only, pair)	\$425
<input type="checkbox"/>	PY5660	Elbow Stop with Arm Rest-9"-14.5" Range (for planar and contoured backs, pair)	\$384

Vest and Straps

<input type="checkbox"/>	PY5620	Chest Strap (circumference up to 42", velcro with D-ring)	\$159
<input type="checkbox"/>	PY5634	High Mount Chest Vest Bracket (head support required)	\$84
<input type="checkbox"/>	PNG30170	X-Style Chest Vest-11"Lx9.5"W	\$297
<input type="checkbox"/>	PNG30056	X-Style Chest Vest-14.5"Lx10.5"W	\$297
<input type="checkbox"/>	PNG30171	Y-Style Chest Vest-17.5"Lx11.5"W	\$297
<input type="checkbox"/>	PNG30021	Y-Style Chest Vest-20.5"Lx12.5"W	\$297
<input type="checkbox"/>	PNG30318	Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$345
<input type="checkbox"/>	PNG30322	Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W	\$345
<input type="checkbox"/>	PNG30320	Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W	\$345
<input type="checkbox"/>	PNG30324	Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W	\$345

Head Supports-supine requires a head support

<input type="checkbox"/>	PY5626	Head Support-5"Hx8"W	\$388
<input type="checkbox"/>	PY5628	Head Support-6"Hx10"W	\$388
<input type="checkbox"/>	PY5630	Head Support-7"Hx14"W	\$388
<input type="checkbox"/>	PA5586	Head Support-Form to Fit-5"Hx10"W	\$388

Other Accessories

<input type="checkbox"/>	P80809	Tool Pouch	\$27
--------------------------	--------	------------	------

Notes:

(A) Seat Depth Range: 15"-20"
(seat pivot to planar back)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 7"-18.5"
(from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

info@easystand.com easystand.com

Email form to webquotes@easystand.com