

**Bill To** (Ultimate Medical Supplier Only)

Purchase Order #:	
Quote #:	
Account #:	
Ordered By:	
Supplier Company:	
Address:	
City, State, Zip:	
Phone:	Email:
Payer Source:	
Name of ATP (Required):	

Ship To

Name:	
Address:	
City, State, Zip:	
Phone:	
Tag For:	

Bantam Extra Small User Specifications

Approximate Height Range	28"-40"
Weight Limit	50 lbs.
Seat Depth Range (seat to back surface)	7"-12"
Seat to Foot Plate Range (from seat pivot)	4.5"-16.5"

Base *Note: The DEFAULT items will be added if no selection is made.***■ PT50001-1 EasyStand® Bantam Extra Small Base** \$3,138.67**Lift Mechanisms-must select one bolded option below**

- | | | |
|---|----------------|------------|
| <input type="checkbox"/> PT50126 Gas Spring Lift | DEFAULT | N/C |
| <input type="checkbox"/> PT50232 Non-Supine | | N/C |
| <input type="checkbox"/> PT50234 Supine (head support required)* | | \$723.36 |
| SUPINE CANNOT BE ORDERED OFF-UNIT AND ONLY AVAILABLE WITH PT50126 GAS SPRING LIFT | | |
| <input type="checkbox"/> PA5614 Gas Spring Lift Lockout | | \$182.21 |
| <input type="checkbox"/> PT50236 Dual Control (foot and hand control. Push handle required) | | \$182.21 |
| <input type="checkbox"/> PT50020 Manual Hydraulic Actuator with Removable Handle | | \$1,959.10 |
| <input type="checkbox"/> PNG50044 Quad Grip Handle Extension | | \$397.30 |
| <input type="checkbox"/> PNG50070 T-Style Handle Extension | | \$247.97 |
| <input type="checkbox"/> PT50022 Pow'r Up Lift (battery powered) | | \$2,022.12 |

Frame Color-must select one

- | | |
|---|-----|
| <input type="checkbox"/> Green DEFAULT | N/C |
| <input type="checkbox"/> Purple | N/C |
| <input type="checkbox"/> White | N/C |
| <input type="checkbox"/> Orange | N/C |
- Go to easystand.com to see specific frame colors*

Front Frame-must select one

- | | | |
|--|----------------|------------|
| <input type="checkbox"/> PT50016-1 Front Frame (5" Front Wheels) | DEFAULT | \$280.85 |
| <input type="checkbox"/> PT50116 Front Frame (3" Front Swivel Casters) | | \$350.72 |
| <input type="checkbox"/> PT50177 Mobile (not available with supine) | | \$3,090.72 |

Trays/Anterior Support-must select one

- | | |
|---|------------|
| <input type="checkbox"/> PA5602 Clear Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable)* | \$1,412.47 |
| <input type="checkbox"/> PA5600 Black Molded Swing-Away Tray-11"Lx19"W DEFAULT* | \$1,190.53 |
| (depth, height & angle adjustable) | |
| <input type="checkbox"/> PT50204 Black Molded Tray-11"Lx19"W (depth & height adjustable)* | \$775.42 |
| <input type="checkbox"/> PA3038 Hygienic Padded Tray Cover (11"Lx19"W tray only, must select pattern) | \$347.98 |
| <input type="checkbox"/> PT50076 Clear Tray-21"Lx24"W (depth, height & angle adjustable)* | \$1,398.77 |
| <input type="checkbox"/> PT50074 Black Molded Tray-21"Lx24"W (depth, height & angle adjustable)* | \$968.59 |
| <input type="checkbox"/> PT30093 Hygienic Padded Tray Cover(21"Lx24"W tray only, must select pattern) | \$402.78 |
| <input type="checkbox"/> PT50246 Mobile Tray-4"Lx12"W | \$786.38 |
| (depth and height adjustable, includes mobile chest pad, not available w/form to fit backs, lateral supports or elbow stop) | |

Chest Pad*-must select one

- | | |
|---|----------|
| <input type="checkbox"/> PA3123 Pivoting Chest Pad | \$112.34 |
| <input type="checkbox"/> PT30046 Small Chest Pad DEFAULT | N/C |

Tray Accessories*

- | | |
|---|----------|
| <input type="checkbox"/> PNG50042 Hand Grips (for clear tray only, pair) | \$304.14 |
| <input type="checkbox"/> PT50156 Elbow Pad-9.25" Cutout (21"Lx24"W tray only) | \$241.12 |
| <input type="checkbox"/> PT50160 Elbow Pad-11" Cutout (21"Lx24"W tray only) | \$257.56 |
| <input type="checkbox"/> PA5604 Forearm Wings-Black Molded Tray (11"Lx19"W tray only, pair) | \$241.12 |
| <input type="checkbox"/> PA5606 Forearm Wings-Clear Tray (11"Lx19"W tray only, pair) | \$241.12 |

Knee Supports-must select one

- | | |
|---|----------|
| <input type="checkbox"/> PT50262-1 Slide-In Knee Pads (pair) DEFAULT | \$528.82 |
| <input type="checkbox"/> PT50264-1 Swing-Away Knee Pads (pair) | \$874.06 |

Knee Pad Size-must select one, inside width

- | | |
|---|----------|
| <input type="checkbox"/> PT50266 2.5" | N/C |
| <input type="checkbox"/> PT50268 3.25" DEFAULT | N/C |
| <input type="checkbox"/> PT50270 4.25" | N/C |
| <input type="checkbox"/> PT50322 Asymmetric-2.5" | \$57.54 |
| <input type="checkbox"/> PT50324 Asymmetric-3.25" | \$68.50 |
| <input type="checkbox"/> PT50326 Asymmetric-4.25" | \$112.34 |

Knee Pad Accessories

- | | |
|---|----------|
| <input type="checkbox"/> PT50278 Calf Pad | \$241.12 |
|---|----------|

Foot Plates-must select one

- | | |
|---|----------|
| <input type="checkbox"/> PT50250-1 Height Adjustable Foot Plates (pair) | \$427.44 |
| <input type="checkbox"/> PT50252 Multi-Adjustable Foot Plates (pair) DEFAULT | \$531.56 |

Foot Holder Size-must select one

- | | |
|---|-----|
| <input type="checkbox"/> PT50254 6"Lx3"W | N/C |
| <input type="checkbox"/> PT50256 7.75"Lx3.25"W DEFAULT | N/C |
| <input type="checkbox"/> PT50258 9.75"Lx4"W | N/C |

* not available with mobile



Foot Plate Accessories

<input type="checkbox"/> PT50080	Foot Straps-10"L (length over top of foot, pair)	\$110.97
<input type="checkbox"/> PY5636	Foot Straps-15"L (length over top of foot, pair)	\$126.04
<input type="checkbox"/> PT50082	Secure Foot Straps-10"L (length over top of foot, two pair)	\$189.06
<input type="checkbox"/> PY5638	Secure Foot Straps-15"L (length over top of foot, two pair)	\$202.76

Seats-must select one

<input type="checkbox"/> PT50042-1	Planar Seat DEFAULT	\$458.95
<input type="checkbox"/> PT30082	Hygienic Cover for Planar Seat (must select pattern)	\$130.15
<input type="checkbox"/> PT50044	Form to Fit Seat	\$569.92
<input type="checkbox"/> PT30083	Hygienic Cover for Form to Fit Seat (must select pattern)	\$130.15

Seat Accessories

<input type="checkbox"/> PNG50499	Easy-Adjust Seat Depth- CANNOT be ordered off-unit (roller bearing allows for smoother seat depth adjustment)	\$175.36
-----------------------------------	---	----------

Pelvic Support

<input type="checkbox"/> PT50098	Hip Supports (4"-8" range, pair, not available with form to fit seat)	\$583.62
<input type="checkbox"/> PA5840	Hygienic Cover for Hip Supports (must select pattern, pair)	\$168.51
<input type="checkbox"/> PT30066	Positioning Belt (fits hip circumference up to 23", velcro with D-ring)	\$110.97
<input type="checkbox"/> PT30094	Positioning Belt (fits hip circumference up to 33", velcro with D-ring)	\$119.19

Backs-height from seat, must select one

<input type="checkbox"/> PT50050-1	Planar Back-11"-13"H Range DEFAULT	\$458.95
<input type="checkbox"/> PT30078	Hygienic Cover for PT50050-1 (must select pattern)	\$130.15
<input type="checkbox"/> PT50238-1	Planar Back-15"-17"H Range	\$480.87
<input type="checkbox"/> PT30089	Hygienic Cover for PT50238-1 (must select pattern)	\$130.15
<input type="checkbox"/> PT50052	Form to Fit Back-11"-13"H Range	\$569.92
<input type="checkbox"/> PT30079	Hygienic Cover for PT50052 (must select pattern)	\$130.15
<input type="checkbox"/> PT50240	Form to Fit Back-15"-17"H Range	\$593.21
<input type="checkbox"/> PT30090	Hygienic Cover for PT50240 (must select pattern)	\$130.15

Lateral Supports*

<input type="checkbox"/> PT50332	Lateral Supports-5"-9" Range Square Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PT50330	Lateral Supports-8"-12" Range Square Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PA5840	Hygienic Cover for Square Pads (must select pattern, pair)	\$168.51
<input type="checkbox"/> PT50248	Lateral Supports-5"-9" Range Curved Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PT50008	Lateral Supports-8"-12" Range Curved Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PT50276	Lateral Supports-5"-9" Range Flat Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PT50274	Lateral Supports-8"-12" Range Flat Pad (for planar back only, pair)	\$583.62
<input type="checkbox"/> PB5720	Hygienic Cover for Curved/Flat Pad (must select pattern, pair)	\$171.25
<input type="checkbox"/> PT50216	Lateral Supports with Elbow Stop-8"-12" Range (for planar back only, pair)	\$597.32
<input type="checkbox"/> PT50218	Lateral Supports with Elbow Stop and Arm Rest-8"-12" Range (for planar back only, pair)	\$660.34
<input type="checkbox"/> PT50220	Elbow Stop with Arm Rest-8"-12" Range (for planar and form to fit back, pair)	\$597.32

Upper Body Support

<input type="checkbox"/> PA5672	Padded Positioning Strap	\$247.97
<input type="checkbox"/> PT50086	X-Style Chest Vest-9.5"Lx9"W	\$461.69
<input type="checkbox"/> PT50088	X-Style Chest Vest-11"Lx9.5"W	\$476.76
<input type="checkbox"/> PT50302	Non-Stretch X-Style Chest Vest-9.5"Lx9"W	\$537.04
<input type="checkbox"/> PNG30318	Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$550.74
<input type="checkbox"/> PT50128	High Mount Chest Vest Bracket	\$127.41

Head Support-supine requires a head support

<input type="checkbox"/> PT50062	Head Support-5"Hx8"W	\$602.80
<input type="checkbox"/> PA3035	Hygienic Cover for Head Support-5"Hx8"W (must select pattern)	\$168.51
<input type="checkbox"/> PT50064	Head Support-6"Hx10"W	\$617.87
<input type="checkbox"/> PA3036	Hygienic Cover for Head Support-6"Hx10"W (must select pattern)	\$168.51
<input type="checkbox"/> PT50284	Head Support-Form to Fit-5"Hx10"W	\$649.38
<input type="checkbox"/> PA3037	Hygienic Cover for Head Support-Form to Fit (must select pattern)	\$168.51

Other Accessories

<input type="checkbox"/> P80809	Tool Pouch	\$45.21
<input type="checkbox"/> PT50212	Push Handle for Supine (required for dual control)	\$247.97
<input type="checkbox"/> PT50230	Push Handle for Non-Supine (required for dual control)	\$247.97
<input type="checkbox"/> P80246	Angle Locator	\$63.02

Hygienic Cover Patterns (must select one if hygienic covers were ordered)

- ☐ Rainbow Unicorn
- ☐ Sea Creatures
- ☐ Fancy Pants Elephants
- ☐ Space Race
- ☐ All Smiles Animals
- ☐ Black (not available on hip support pads, lateral pads or padded tray cover)

Go to easystand.com to see specific fabric patterns

* not available with mobile

Notes:

(A) Seat Depth Range: 7"-12" (seat to back surface)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 4.5"-16.5" (from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



AltimateMedical
creating products · changing lives

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

webquotes@altimatemedical.com easystand.com



Email form to webquotes@altimatemedical.com

