

**Bill To** (Ultimate Medical Supplier Only)

Purchase Order #: _____
Quote #: _____
Account #: _____
Ordered By: _____
Supplier Company: _____
Address: _____
City, State, Zip: _____
Phone: _____ Email: _____
Payer Source: _____
Name of ATP (Required): _____

Ship To

Name: _____
Address: _____
City, State, Zip: _____
Phone: _____
Tag For: _____

Bantam Extra Small User Specifications

Approximate Height Range 28"-40"
Weight Limit 50 lbs.
Seat Depth Range (seat to back surface) 7"-12"
Seat to Foot Plate Range (from seat pivot) 4.5"-16.5"

Base *Note: The DEFAULT items will be added if no selection is made.*

☐ **PT50001-1 EasyStand® Bantam Extra Small Base** \$2,291

Lift Mechanisms-must select one bolded option below

☐ **PT50126 Gas Spring Lift** **DEFAULT** N/C
☐ PT50232 Non-Supine N/C
☐ PT50234 Supine (head support required)* \$528

SUPINE CANNOT BE ORDERED OFF-UNIT AND ONLY AVAILABLE WITH PT50126 GAS SPRING LIFT

☐ PA5614 Gas Spring Lift Lockout \$133
☐ PT50236 Dual Control (foot and hand control. Push handle required) \$133

☐ **PT50020 Manual Hydraulic Actuator with Removable Handle** \$1,430
☐ PNG50044 Quad Grip Handle Extension \$290
☐ PNG50070 T-Style Handle Extension \$181

☐ **PT50022 Pow'r Up Lift** (battery powered) \$1,476

Frame Color-must select one

☐ Green **DEFAULT** N/C
☐ Purple N/C
☐ White N/C
☐ Orange N/C

Go to easystand.com to see specific frame colors

Front Frame-must select one

☐ PT50016-1 Front Frame (5" Front Wheels) **DEFAULT** \$205
☐ PT50116 Front Frame (3" Front Swivel Casters) \$256
☐ PT50177 Mobile (not available with supine) \$2,256

Trays/Anterior Support-must select one

☐ PA5602 Clear Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable)* \$1,031
☐ PA5600 Black Molded Swing-Away Tray-11"Lx19"W **DEFAULT*** \$869
(depth, height & angle adjustable)
☐ PT50204 Black Molded Tray-11"Lx19"W (depth & height adjustable)* \$566
☐ PA3038 Hygienic Padded Tray Cover (11"Lx19"W tray only, must select pattern) \$254
☐ PT50076 Clear Tray-21"Lx24"W (depth, height & angle adjustable)* \$1,021
☐ PT50074 Black Molded Tray-21"Lx24"W (depth, height & angle adjustable)* \$707
☐ PT30093 Hygienic Padded Tray Cover (21"Lx24"W tray only, must select pattern) \$294
☐ PT50246 Mobile Tray-4"Lx12"W \$574
(depth and height adjustable, includes mobile chest pad, not available w/form to fit backs, lateral supports or elbow stop)

Chest Pad*-must select one

☐ PA3123 Pivoting Chest Pad \$82
☐ PT30046 Small Chest Pad **DEFAULT** N/C

Tray Accessories*

☐ PNG50042 Hand Grips (for clear tray only, pair) \$222
☐ PT50156 Elbow Pad-9.25" Cutout (21"Lx24"W tray only) \$176
☐ PT50160 Elbow Pad-11" Cutout (21"Lx24"W tray only) \$188
☐ PA5604 Forearm Wings-Black Molded Tray (11"Lx19"W tray only, pair) \$176
☐ PA5606 Forearm Wings-Clear Tray (11"Lx19"W tray only, pair) \$176

Knee Supports-must select one

☐ PT50262-1 Slide-In Knee Pads (pair) **DEFAULT** \$386
☐ PT50264-1 Swing-Away Knee Pads (pair) \$638

Knee Pad Size-must select one, inside width

☐ PT50266 2.5" N/C
☐ PT50268 3.25" **DEFAULT** N/C
☐ PT50270 4.25" N/C
☐ PT50322 Asymmetric-2.5" \$42
☐ PT50324 Asymmetric-3.25" \$50
☐ PT50326 Asymmetric-4.25" \$82

Knee Pad Accessories

☐ PT50278 Calf Pad \$176

Foot Plates-must select one

☐ PT50250-1 Height Adjustable Foot Plates (pair) \$312
☐ PT50252 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$388

Foot Holder Size-must select one

☐ PT50254 6"Lx3"W N/C
☐ PT50256 7.75"Lx3.25"W **DEFAULT** N/C
☐ PT50258 9.75"Lx4"W N/C

* not available with mobile

Foot Plate Accessories

<input type="checkbox"/> PT50080	Foot Straps-10"L (length over top of foot, pair)	\$81
<input type="checkbox"/> PY5636	Foot Straps-15"L (length over top of foot, pair)	\$92
<input type="checkbox"/> PT50082	Secure Foot Straps-10"L (length over top of foot, two pair)	\$138
<input type="checkbox"/> PY5638	Secure Foot Straps-15"L (length over top of foot, two pair)	\$148

Seats-must select one

<input type="checkbox"/> PT50042-1	Planar Seat DEFAULT	\$335
<input type="checkbox"/> PT30082	Hygienic Cover for Planar Seat (must select pattern)	\$95
<input type="checkbox"/> PT50044	Form to Fit Seat	\$416
<input type="checkbox"/> PT30083	Hygienic Cover for Form to Fit Seat (must select pattern)	\$95

Seat Accessories

<input type="checkbox"/> PNG50499	Easy-Adjust Seat Depth- CANNOT be ordered off-unit (roller bearing allows for smoother seat depth adjustment)	\$128
-----------------------------------	-------------------------------------------------------------------------------------------------------------------------	-------

Pelvic Support

<input type="checkbox"/> PT50098	Hip Supports (4"-8" range, pair, not available with form to fit seat)	\$426
<input type="checkbox"/> PA5840	Hygienic Cover for Hip Supports (must select pattern, pair)	\$123
<input type="checkbox"/> PT30066	Positioning Belt (fits hip circumference up to 23", velcro with D-ring)	\$81
<input type="checkbox"/> PT30094	Positioning Belt (fits hip circumference up to 33", velcro with D-ring)	\$87

Backs-height from seat, must select one

<input type="checkbox"/> PT50050-1	Planar Back-11"-13"H Range DEFAULT	\$335
<input type="checkbox"/> PT30078	Hygienic Cover for PT50050-1 (must select pattern)	\$95
<input type="checkbox"/> PT50238-1	Planar Back-15"-17"H Range	\$351
<input type="checkbox"/> PT30089	Hygienic Cover for PT50238-1 (must select pattern)	\$95
<input type="checkbox"/> PT50052	Form to Fit Back-11"-13"H Range	\$416
<input type="checkbox"/> PT30079	Hygienic Cover for PT50052 (must select pattern)	\$95
<input type="checkbox"/> PT50240	Form to Fit Back-15"-17"H Range	\$433
<input type="checkbox"/> PT30090	Hygienic Cover for PT50240 (must select pattern)	\$95

Lateral Supports*

<input type="checkbox"/> PT50332	Lateral Supports-5"-9" Range Square Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PT50330	Lateral Supports-8"-12" Range Square Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PA5840	Hygienic Cover for Square Pads (must select pattern, pair)	\$123
<input type="checkbox"/> PT50248	Lateral Supports-5"-9" Range Curved Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PT50008	Lateral Supports-8"-12" Range Curved Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PT50276	Lateral Supports-5"-9" Range Flat Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PT50274	Lateral Supports-8"-12" Range Flat Pad (for planar back only, pair)	\$426
<input type="checkbox"/> PB5720	Hygienic Cover for Curved/Flat Pad (must select pattern, pair)	\$125
<input type="checkbox"/> PT50216	Lateral Supports with Elbow Stop-8"-12" Range (for planar back only, pair)	\$436
<input type="checkbox"/> PT50218	Lateral Supports with Elbow Stop and Arm Rest-8"-12" Range (for planar back only, pair)	\$482
<input type="checkbox"/> PT50220	Elbow Stop with Arm Rest-8"-12" Range (for planar and form to fit back, pair)	\$436

Upper Body Support

<input type="checkbox"/> PA5672	Padded Positioning Strap	\$181
<input type="checkbox"/> PT50086	X-Style Chest Vest-9.5"Lx9"W	\$337
<input type="checkbox"/> PT50088	X-Style Chest Vest-11"Lx9.5"W	\$348
<input type="checkbox"/> PT50302	Non-Stretch X-Style Chest Vest-9.5"Lx9"W	\$392
<input type="checkbox"/> PNG30318	Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$402
<input type="checkbox"/> PT50128	High Mount Chest Vest Bracket	\$93

Head Support-supine requires a head support

<input type="checkbox"/> PT50062	Head Support-5"Hx8"W	\$440
<input type="checkbox"/> PA3035	Hygienic Cover for Head Support-5"Hx8"W (must select pattern)	\$123
<input type="checkbox"/> PT50064	Head Support-6"Hx10"W	\$451
<input type="checkbox"/> PA3036	Hygienic Cover for Head Support-6"Hx10"W (must select pattern)	\$123
<input type="checkbox"/> PT50284	Head Support-Form to Fit-5"Hx10"W	\$474
<input type="checkbox"/> PA3037	Hygienic Cover for Head Support-Form to Fit (must select pattern)	\$123

Other Accessories

<input type="checkbox"/> P80809	Tool Pouch	\$33
<input type="checkbox"/> PT50212	Push Handle for Supine (required for dual control)	\$181
<input type="checkbox"/> PT50230	Push Handle for Non-Supine (required for dual control)	\$181
<input type="checkbox"/> P80246	Angle Locator	\$46

Hygienic Cover Patterns (must select one if hygienic covers were ordered)

<input type="checkbox"/> Rainbow Unicorn	
<input type="checkbox"/> Sea Creatures	<i>Go to easystand.com to see specific fabric patterns</i>
<input type="checkbox"/> Fancy Pants Elephants	
<input type="checkbox"/> Space Race	
<input type="checkbox"/> All Smiles Animals	
<input type="checkbox"/> Black (not available on hip support pads, lateral pads or padded tray cover)	

* not available with mobile

Notes:

(A) Seat Depth Range: 7"-12"
(seat to back surface)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 4.5"-16.5"
(from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



AltimateMedical
creating products · changing lives

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

webquotes@altimatemedical.com easystand.com

Email form to webquotes@altimatemedical.com