



Bill To (Ultimate Medical Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____ Email: _____

Payer Source: _____

Name of ATP (Required): _____

Ship To

Name: _____

Address: _____

City, State, Zip: _____

Phone: _____

Tag For: _____

Bantam Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat pivot to planar back)	15"-20"
Seat to Foot Plate Range (from seat pivot)	7"-18.5"

Base-must select one *Note: The DEFAULT items will be added if no selection is made.*

- PY5690 EasyStand® Bantam Medium Supine Ready Base** \$2,522
 - Supine-must select one**
 - Supine CANNOT be ordered off-unit, Hand Control right side only**
 - PY5688 Supine (head support required) **DEFAULT** \$758
 - PY5689 Non-Supine N/C

- PY5691 EasyStand® Bantam Medium Base** \$2,522
 - Supine not available with PY5691**
 - Choose actuator handle location-only for PY5691**
 - PY5698 Right Side **DEFAULT** N/C
 - PY5699 Left Side N/C

- Lift Mechanisms Accessories**
- PY5632 Adjustable/Removable Actuator Handle (not available with PY5690) \$268
 - PNG50044 Quad Grip Handle Extension (not available with PY5690) \$290
 - PNG50070 T-Style Handle Extension (not available with PY5690) \$181

- Front Frame-must select one**
- PY5694 Front Frame (5" Front Wheels) **DEFAULT** \$244
 - PY5695 Front Frame (3" Front Swivel Casters) \$301

- Shadow Trays-must select one**
- PY5604 Black Molded Tray-21"Lx24"W (depth, height and angle adjustable) \$776
 - PY5606 Clear Tray-21"Lx24"W (depth, height and angle adjustable) \$1,073
 - PY5656 Black Molded Swing-Away Tray-21"Lx24"W **DEFAULT** (depth, height and angle adjustable) \$1,034
 - PY5657 Clear Swing-Away Tray-21"Lx24"W (depth, height and angle adjustable) \$1,256
 - PNG30295 Padded Tray Cover (not available with elbow pads) \$266

- Chest Pad-must select one**
- PA3123 Pivoting Chest Pad \$82
 - PT30046 Small Chest Pad **DEFAULT** N/C
 - PNG30000 Large Contoured Chest Pad (not available with elbow pads) \$168

- Tray Accessories**
- PNG50042 Hand Grips (for clear trays only, pair) \$222
 - PT50156 Elbow Pad-9.25" Cutout \$176
 - PT50160 Elbow Pad-11" Cutout \$188
 - PY5582 Elbow Pad-13" Cutout \$199

- Knee Pads-must select one**
- PY5594-1 Slide-In Knee Pads (pair) \$493
 - PY5596-1 Swing-Away Knee Pads (pair) **DEFAULT** \$794

- Knee Pad Size-must select one**
- PY5624 4.0" N/C
 - PY5598 5.5" **DEFAULT** N/C
 - PY5600-1 5.5" Roho® Knee Pads (pair) \$510

- Foot Plates-must select one**
- PY5584-1 Height-Adjustable Foot Plates (pair) \$349
 - PY5586-1 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$388

- Foot Plate Size-must select one**
- PY5588 7.75"Lx3.25"W (pair) N/C
 - PY5590 9.75"Lx4"W (pair) N/C
 - PY5592 11.75"Lx5"W (pair) **DEFAULT** N/C

- Foot Plate Accessories**
- PT50080 Foot Straps-10"L (length over top of foot, pair) \$81
 - PY5636 Foot Straps-15"L (length over top of foot, pair) \$92
 - PT50082 Secure Foot Straps-10"L (length over top of foot, two pair) \$138
 - PY5638 Secure Foot Straps-15"L (length over top of foot, two pair) \$148

Seats-must select one

<input type="checkbox"/> PY5522-1	Planar Seat DEFAULT	\$527
<input type="checkbox"/> PY3042	Black Hygienic Seat Cover for PY5522-1	\$145
<input type="checkbox"/> PY5570	Contoured Seat-Narrow 11"-13"W	\$510
<input type="checkbox"/> PY3043	Black Hygienic Seat Cover for PY5570	\$145
<input type="checkbox"/> PY5572	Contoured Seat-Wide 12"-14"W	\$527
<input type="checkbox"/> PY3044	Black Hygienic Seat Cover for PY5572	\$145

Seat Accessories

<input type="checkbox"/> PY5645	Easy-Adjust Seat Depth- CANNOT be ordered off-unit (roller bearing allows for smoother seat depth adjustment)	\$128
---------------------------------	---	-------

Pelvic Support

<input type="checkbox"/> PY5560	Hip Supports (9"-15" range, pair, not available with contoured seat)	\$426
<input type="checkbox"/> PY3008	Velcro® Positioning Belt (fits hip circumference up to 66", velcro with D-ring)	\$97
<input type="checkbox"/> PY3014	Positioning Belt with Airline Style Buckle (fits hip circumference up to 66")	\$199

Backs-must select one

<input type="checkbox"/> PY5530	Planar Back-16"H	\$857
<input type="checkbox"/> PY3041	Black Hygienic Back Cover for PY5530	\$108
<input type="checkbox"/> PY5532	Planar Back-20"H DEFAULT	\$905
<input type="checkbox"/> PY3040	Black Hygienic Back Cover for PY5532	\$108
<input type="checkbox"/> PY5574	Contoured Back-16"H	\$945
<input type="checkbox"/> PY3045	Black Hygienic Back Cover for PY5574	\$108
<input type="checkbox"/> PY5576	Contoured Back-20"H	\$993
<input type="checkbox"/> PY3046	Black Hygienic Back Cover for PY5576	\$108

Back Accessories

<input type="checkbox"/> PNG50040	Push Handles (pair)	\$181
<input type="checkbox"/> PY5643	Angle Locator	\$46

Lateral Supports

<input type="checkbox"/> PY5686	Lateral Supports with Flat Pads-9"-14.5" Range (for planar back only, pair)	\$426
<input type="checkbox"/> PY5562	Lateral Supports with Curved Pads-9"-14.5" Range (for planar back only, pair)	\$426
<input type="checkbox"/> PY5658	Lateral Supports with Elbow Stop-9"-14.5" Range (for planar back only, pair)	\$436
<input type="checkbox"/> PY5659	Lateral Supports with Elbow Stop and Arm Rest-9"-14.5" Range (for planar back only, pair)	\$510
<input type="checkbox"/> PY5660	Elbow Stop with Arm Rest-9"-14.5" Range (for planar and contoured backs, pair)	\$436

Vest and Straps

<input type="checkbox"/> PY5620	Chest Strap (circumference up to 42", velcro with D-ring)	\$204
<input type="checkbox"/> PY5634	High Mount Chest Vest Bracket (head support required)	\$96
<input type="checkbox"/> PNG30170	X-Style Chest Vest-11"Lx9.5"W	\$348
<input type="checkbox"/> PNG30056	X-Style Chest Vest-14.5"Lx10.5"W	\$355
<input type="checkbox"/> PNG30171	Y-Style Chest Vest-17.5"Lx11.5"W	\$348
<input type="checkbox"/> PNG30021	Y-Style Chest Vest-20.5"Lx12.5"W	\$355
<input type="checkbox"/> PNG30318	Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$402
<input type="checkbox"/> PNG30322	Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W	\$409
<input type="checkbox"/> PNG30320	Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W	\$402
<input type="checkbox"/> PNG30324	Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W	\$409

Head Supports-supine requires a head support

<input type="checkbox"/> PY5626	Head Support-5"Hx8"W	\$440
<input type="checkbox"/> PY5628	Head Support-6"Hx10"W	\$451
<input type="checkbox"/> PY5630	Head Support-7"Hx14"W	\$463
<input type="checkbox"/> PA5586	Head Support-Form to Fit-5"Hx10"W	\$474

Other Accessories

<input type="checkbox"/> P80809	Tool Pouch	\$33
---------------------------------	------------	------

Notes:

(A) Seat Depth Range: 15"-20" (seat pivot to planar back)

Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 7"-18.5" (from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



AltimateMedical
creating products · changing lives

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

webquotes@altimatemedical.com easystand.com

Email form to webquotes@altimatemedical.com