

| Bill To (Altimate Medical Supplier Only) |        |
|------------------------------------------|--------|
| Purchase Order #:                        |        |
| Quote #:                                 |        |
| Account #:                               |        |
| Ordered By:                              |        |
| Supplier Company:                        |        |
| Address:                                 |        |
| City, State, Zip:                        |        |
| Phone:                                   | Email: |
| Payer Source:                            |        |
| Name of ATP (Required):                  |        |
| Ship To                                  |        |
| Name:                                    |        |
| Address:                                 |        |
| City, State, Zip:                        |        |
| Phone:                                   |        |
| Tag For:                                 |        |

## Bantam Small User Specifications

|                                            |            |
|--------------------------------------------|------------|
| Approximate Height Range                   | 36"-54"    |
| Weight Limit                               | 100 lbs.   |
| Seat Depth Range (seat to back surface)    | 11"-16"    |
| Seat to Foot Plate Range (from seat pivot) | 4.5"-16.5" |

**Base** *Note: The DEFAULT items will be added if no selection is made.*

☒ **PT50002-1 EasyStand® Bantam Small Base** \$2,291

## Lift Mechanisms-must select one bolded option below

- ☐ **PT50018 Gas Spring Lift** **DEFAULT** N/C
- ☐ PT50232 Non-Supine N/C
- ☐ PT50234 Supine (head support required)\* \$528
- SUPINE CANNOT BE ORDERED OFF-UNIT AND ONLY AVAILABLE WITH PT50018 GAS SPRING LIFT**
- ☐ PA5614 Gas Spring Lift Lockout \$133
- ☐ PT50236 Dual Control (foot and hand control. Push handle required) \$133
- ☐ **PT50020 Manual Hydraulic Actuator with Removable Handle** \$1,430
- ☐ PNG50044 Quad Grip Handle Extension \$290
- ☐ PNG50070 T-Style Handle Extension \$181
- ☐ **PT50022 Pow'r Up Lift** (battery powered) \$1,476

## Frame Color-must select one

- ☐ Green **DEFAULT** N/C
- ☐ Purple N/C
- ☐ White N/C
- ☐ Orange N/C

*Go to easystand.com to see specific frame colors*

## Front Frame-must select one

- ☐ PT50016-1 Front Frame (5" Front Wheels) **DEFAULT** \$205
- ☐ PT50116 Front Frame (3" Front Swivel Casters) \$256
- ☐ PT50177 Mobile (not available with supine) \$2,256

## Trays/Anterior Support-must select one

- ☐ PA5602 Clear Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable)\* \$1,031
- ☐ PA5600 Black Molded Swing-Away Tray-11"Lx19"W **DEFAULT**\* \$869  
(depth, height & angle adjustable)
- ☐ PT50204 Black Molded Tray-11"Lx19"W (depth & height adjustable)\* \$566
- ☐ PA3038 Hygienic Padded Tray Cover (11"Lx19"W tray only, must select pattern) \$254
- ☐ PT50076 Clear Tray-21"Lx24"W (depth, height & angle adjustable)\* \$1,021
- ☐ PT50074 Black Molded Tray-21"Lx24"W (depth, height & angle adjustable)\* \$707
- ☐ PT30093 Hygienic Padded Tray Cover (21"Lx24"W tray only, must select pattern) \$294
- ☐ PT50246 Mobile Tray-4"Lx12"W \$574  
(depth and height adjustable, includes mobile chest pad, not available w/form to fit backs, lateral supports or elbow stop)

## Chest Pad\*-must select one

- ☐ PA3123 Pivoting Chest Pad \$82
- ☐ PT30046 Small Chest Pad **DEFAULT** N/C

## Tray Accessories\*

- ☐ PNG50042 Hand Grips (for clear tray only, pair) \$222
- ☐ PT50156 Elbow Pad-9.25" Cutout (21"Lx24"W tray only) \$176
- ☐ PT50160 Elbow Pad-11" Cutout (21"Lx24"W tray only) \$188
- ☐ PA5604 Forearm Wings-Black Molded Tray (11"Lx19"W tray only, pair) \$176
- ☐ PA5606 Forearm Wings-Clear Tray (11"Lx19"W tray only, pair) \$176

## Knee Supports-must select one

- ☐ PT50262-1 Slide-In Knee Pads (pair) **DEFAULT** \$386
- ☐ PT50264-1 Swing-Away Knee Pads (pair) \$638

## Knee Pad Size-must select one, inside width

- ☐ PT50268 3.25" N/C
- ☐ PT50270 4.25" **DEFAULT** N/C
- ☐ PT50272 5" N/C
- ☐ PT50324 Asymmetric-3.25" \$50
- ☐ PT50326 Asymmetric-4.25" \$82
- ☐ PT50328 Asymmetric-5" \$88

## Knee Pad Accessories

- ☐ PT50278 Calf Pad \$176

## Foot Plates-must select one

- ☐ PT50250-1 Height Adjustable Foot Plates (pair) \$312
- ☐ PT50252 Multi-Adjustable Foot Plates (pair) **DEFAULT** \$388

## Foot Holder Size-must select one

- ☐ PT50256 7.75"Lx3.25"W **DEFAULT** N/C
- ☐ PT50258 9.75"Lx4"W N/C
- ☐ PT50260 11.75"Lx5"W N/C

\* not available with mobile

## Foot Plate Accessories

|                                  |                                                             |       |
|----------------------------------|-------------------------------------------------------------|-------|
| <input type="checkbox"/> PT50080 | Foot Straps-10"L (length over top of foot, pair)            | \$81  |
| <input type="checkbox"/> PY5636  | Foot Straps-15"L (length over top of foot, pair)            | \$92  |
| <input type="checkbox"/> PT50082 | Secure Foot Straps-10"L (length over top of foot, two pair) | \$138 |
| <input type="checkbox"/> PY5638  | Secure Foot Straps-15"L (length over top of foot, two pair) | \$148 |

## Seats-must select one

|                                    |                                                           |       |
|------------------------------------|-----------------------------------------------------------|-------|
| <input type="checkbox"/> PT50046-1 | Planar Seat <b>DEFAULT</b>                                | \$363 |
| <input type="checkbox"/> PT30084   | Hygienic Cover for Planar Seat (must select pattern)      | \$95  |
| <input type="checkbox"/> PT50048   | Form to Fit Seat                                          | \$445 |
| <input type="checkbox"/> PT30085   | Hygienic Cover for Form to Fit Seat (must select pattern) | \$95  |

## Seat Accessories

|                                   |                                                                                                                         |       |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> PNG50499 | Easy-Adjust Seat Depth- <b>CANNOT be ordered off-unit</b><br>(roller bearing allows for smoother seat depth adjustment) | \$128 |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------|-------|

## Pelvic Support

|                                  |                                                                         |       |
|----------------------------------|-------------------------------------------------------------------------|-------|
| <input type="checkbox"/> PT50007 | Hip Supports (8"-12" range, pair, not available with form to fit seat)  | \$426 |
| <input type="checkbox"/> PA5840  | Hygienic Cover for Hip Supports (must select pattern, pair)             | \$123 |
| <input type="checkbox"/> PT30066 | Positioning Belt (fits hip circumference up to 29", velcro with D-ring) | \$81  |
| <input type="checkbox"/> PT30094 | Positioning Belt (fits hip circumference up to 39", velcro with D-ring) | \$87  |

## Backs-height from seat, must select one

|                                    |                                           |       |
|------------------------------------|-------------------------------------------|-------|
| <input type="checkbox"/> PT50054-1 | Planar Back-14"-16"H Range <b>DEFAULT</b> | \$363 |
| <input type="checkbox"/> PT30080   | Hygienic Cover for PT50054-1              | \$95  |
| <input type="checkbox"/> PT50242-1 | Planar Back-18"-20"H Range                | \$381 |
| <input type="checkbox"/> PT30091   | Hygienic Cover for PT50242-1              | \$95  |
| <input type="checkbox"/> PT50056   | Form to Fit Back-12.5"-14.5"H Range       | \$445 |
| <input type="checkbox"/> PT30081   | Hygienic Cover for PT50056                | \$95  |
| <input type="checkbox"/> PT50244   | Form to Fit Back-16.5"-18.5"H Range       | \$473 |
| <input type="checkbox"/> PT30092   | Hygienic Cover for PT50244                | \$95  |

## Lateral Supports\*

|                                  |                                                                                            |       |
|----------------------------------|--------------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> PT50332 | Lateral Supports-5"-9" Range Square Pad (for planar back only, pair)                       | \$426 |
| <input type="checkbox"/> PT50330 | Lateral Supports-8"-12" Range Square Pad (for planar back only, pair)                      | \$426 |
| <input type="checkbox"/> PA5840  | Hygienic Cover for Square Pads (must select pattern, pair)                                 | \$123 |
| <input type="checkbox"/> PT50248 | Lateral Supports-5"-9" Range Curved Pad (for planar back only, pair)                       | \$426 |
| <input type="checkbox"/> PT50008 | Lateral Supports-8"-12" Range Curved Pad (for planar back only, pair)                      | \$426 |
| <input type="checkbox"/> PT50276 | Lateral Supports-5"-9" Range Flat Pad (for planar back only, pair)                         | \$426 |
| <input type="checkbox"/> PT50274 | Lateral Supports-8"-12" Range Flat Pad (for planar back only, pair)                        | \$426 |
| <input type="checkbox"/> PB5720  | Hygienic Cover for Curved/Flat Pad (must select pattern, pair)                             | \$125 |
| <input type="checkbox"/> PT50216 | Lateral Supports with Elbow Stop-8"-12" Range (for planar back only, pair)                 | \$436 |
| <input type="checkbox"/> PT50218 | Lateral Supports with Elbow Stop and Arm Rest-8"-12" Range<br>(for planar back only, pair) | \$482 |
| <input type="checkbox"/> PT50220 | Elbow Stop with Arm Rest-8"-12" Range (for planar and form to fit back, pair)              | \$436 |

## Upper Body Support

|                                   |                                           |       |
|-----------------------------------|-------------------------------------------|-------|
| <input type="checkbox"/> PA5672   | Padded Positioning Strap                  | \$181 |
| <input type="checkbox"/> PT50086  | X-Style Chest Vest-9.5"Lx9"W              | \$337 |
| <input type="checkbox"/> PT50088  | X-Style Chest Vest-11"Lx9.5"W             | \$348 |
| <input type="checkbox"/> PT50302  | Non-Stretch X-Style Chest Vest-9.5"Lx9"W  | \$392 |
| <input type="checkbox"/> PNG30318 | Non-Stretch X-Style Chest Vest-11"Lx9.5"W | \$402 |
| <input type="checkbox"/> PT50128  | High Mount Chest Vest Bracket             | \$93  |

## Head Support-supine requires a head support

|                                  |                                                                   |       |
|----------------------------------|-------------------------------------------------------------------|-------|
| <input type="checkbox"/> PT50062 | Head Support-5"Hx8"W                                              | \$440 |
| <input type="checkbox"/> PA3035  | Hygienic Cover for Head Support-5"Hx8"W (must select pattern)     | \$123 |
| <input type="checkbox"/> PT50064 | Head Support-6"Hx10"W                                             | \$451 |
| <input type="checkbox"/> PA3036  | Hygienic Cover for Head Support-6"Hx10"W (must select pattern)    | \$123 |
| <input type="checkbox"/> PT50284 | Head Support-Form to Fit-5"Hx10"W                                 | \$474 |
| <input type="checkbox"/> PA3037  | Hygienic Cover for Head Support-Form to Fit (must select pattern) | \$123 |

## Other Accessories

|                                  |                                                        |       |
|----------------------------------|--------------------------------------------------------|-------|
| <input type="checkbox"/> P80809  | Tool Pouch                                             | \$33  |
| <input type="checkbox"/> PT50212 | Push Handle for Supine (required for dual control)     | \$181 |
| <input type="checkbox"/> PT50230 | Push Handle for Non-Supine (required for dual control) | \$181 |
| <input type="checkbox"/> P80246  | Angle Locator                                          | \$46  |

## Hygienic Cover Patterns (must select one if hygienic covers were ordered)

- ☐ Rainbow Unicorn
- ☐ Sea Creatures
- ☐ Fancy Pants Elephants
- ☐ Space Race
- ☐ All Smiles Animals
- ☐ Black (not available on hip support pads, lateral pads or padded tray cover)

*Go to [easystand.com](http://easystand.com) to see specific fabric patterns*

\* not available with mobile

## Notes:

---



---



---



---

### (A) Seat Depth Range: 11"-16" (seat to back surface)

Measure the distance between the back of the buttocks and the back of the knee.

### (B) Seat to Foot Plate Range: 4.5"-16.5" (from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



**AltimateMedical**  
creating products · changing lives

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

[webquotes@altimatemedical.com](mailto:webquotes@altimatemedical.com) [easystand.com](http://easystand.com)

Email form to [webquotes@altimatemedical.com](mailto:webquotes@altimatemedical.com)