

Retail Order Form

USA Effective January 1, 2026 Date: _____

Evolv Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (seat to back surface)	18"-23"
Seat to Foot Plate Range (from seat pivot)	14"-21"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@altimatemedical.com



Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

webquotes@altimatemedical.com easystand.com

Evolv Large Packages

☐ PK205 Minimum Support Package.....\$6,370

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

☐ Right DEFAULT ☐ Left

Removable Back (select one)

☐ Flat Back-15"H ☐ Flat Back-19"H DEFAULT

☐ Contoured Back-15"H ☐ Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

☐ 26"-55" DEFAULT ☐ 38"-67"

Chest Strap (select one, chest measurement)

☐ 30"-46" ☐ 34"-54" DEFAULT ☐ 46"-66"

☐ PK206 Moderate Support Package.....\$8,782

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

☐ Right DEFAULT ☐ Left

Chest Pad (select one)

☐ Non-Adjustable DEFAULT ☐ Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

☐ 15"L DEFAULT ☐ 21"L

Hip Supports (select one range)

☐ 10"-15" ☐ 12"-17" DEFAULT ☐ 17"-23"

Removable Back (select one)

☐ Flat Back-15"H ☐ Flat Back-19"H DEFAULT

☐ Contoured Back-15"H ☐ Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

☐ 26"-55" DEFAULT ☐ 38"-67"

Chest Strap (select one, chest measurement)

☐ 30"-46" ☐ 34"-54" DEFAULT ☐ 46"-66"

Lateral Supports (select one)

☐ 8"-16.5"W Range (select pad style below)

☐ Flat Pads ☐ Curved Pads

☐ 11"-19.5"W Range (select pad style below)

☐ Flat Pads DEFAULT ☐ Curved Pads

☐ PK207 Maximum Support Package.....\$10,157

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front for Shadow Tray, Adjustable/Removable Actuator Handle, Shadow Tray, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Positioning Belt, Accessories Mounting Bracket, Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Multi-Adjustable Foot Plates, Planar Seat, Front Wheels/Casters, 5" Rear Locking Casters.

Front Wheel/Casters (select one)

☐ 5" Front Wheels DEFAULT ☐ Front Swivel Casters

Actuator Handle Location (select one)

☐ Right DEFAULT ☐ Left

Shadow Tray (select one)

☐ Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)

☐ Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

Chest Pad (select one)

☐ Non-Adjustable DEFAULT ☐ Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

☐ 15"L DEFAULT ☐ 21"L

Hip Supports (select one range)

☐ 10"-15" ☐ 12"-17" DEFAULT ☐ 17"-23"

Removable Back (select one)

☐ Flat Back-15"H ☐ Flat Back-19"H

☐ Contoured Back-15"H ☐ Contoured Back-19"H DEFAULT

Velcro® Positioning Belt (select one, hip circumference)

☐ 26"-55" DEFAULT ☐ 38"-67"

Chest Support (select one)

☐ X-Style 11"Lx9.5"W

☐ X-Style 14.5"Lx10.5"W DEFAULT

☐ Y-Style 17.5"Lx11.5"W

☐ Y-Style 20.5"Lx12.5"W

☐ Non-Stretch X-Style 11"Lx9.5"W

☐ Non-Stretch X-Style 14.5"Lx10.5"W

☐ Non-Stretch Y-Style 17.5"Lx11.5"W ☐ Non-Stretch Y-Style 20.5"Lx12.5"W

Lateral Supports (select one)

☐ 8"-16.5"W Range (select pad style below)

☐ Flat Pads ☐ Curved Pads

☐ 11"-19.5"W Range (select pad style below)

☐ Flat Pads DEFAULT ☐ Curved Pads

Head Support (select one)

☐ Short ☐ Tall DEFAULT

☐ PK208 Minimum Support Package (Mobile)....\$9,289

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

Actuator Handle Location (select one)

☐ Right DEFAULT ☐ Left

Chest Pad (select one)

☐ Non-Adjustable DEFAULT ☐ Adjustable

Removable Back (select one)

☐ Flat Back-15"H ☐ Flat Back-19"H DEFAULT

☐ Contoured Back-15"H ☐ Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

☐ 26"-55" DEFAULT ☐ 38"-67"

Chest Strap (select one, chest measurement)

☐ 30"-46" ☐ 34"-54" DEFAULT ☐ 46"-66"