

## Retail Price List & Order Form

USA Effective May 1, 2021 Date: \_\_\_\_\_

### Evolv Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (seat to back surface)	18"-23"
Seat to Foot Plate Range (from seat pivot)	14"-21"

### Bill To (EasyStand Supplier Only)

Purchase Order #: \_\_\_\_\_

Quote #: \_\_\_\_\_

Account #: \_\_\_\_\_

Ordered By: \_\_\_\_\_

Name of ATP (Required): \_\_\_\_\_

Supplier Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Payer Source: \_\_\_\_\_

### Ship To

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City, State: \_\_\_\_\_

Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Tax For: \_\_\_\_\_

### Email form to [webquotes@easystand.com](mailto:webquotes@easystand.com)

## AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

[easystand.com](http://easystand.com)

# Evolv Large Packages

### PK205 Minimum Support Package.....\$4,307

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

#### Actuator Handle Location (select one)

Right DEFAULT  Left

#### Removable Back (select one)

Flat Back-15"H  Flat Back-19"H DEFAULT

Contoured Back-15"H  Contoured Back-19"H

#### Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT  38"-67"

#### Chest Strap (select one, chest measurement)

30"-46"  34"-54" DEFAULT  46"-66"

### PK206 Moderate Support Package.....\$6,322

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

#### Actuator Handle Location (select one)

Right DEFAULT  Left

#### Chest Pad (select one)

Non-Adjustable DEFAULT  Adjustable

#### Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT  21"L

#### Hip Supports (select one range)

10"-15"  12"-17" DEFAULT  17"-23"

#### Removable Back (select one)

Flat Back-15"H  Flat Back-19"H DEFAULT

Contoured Back-15"H  Contoured Back-19"H

#### Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT  38"-67"

#### Chest Strap (select one, chest measurement)

30"-46"  34"-54" DEFAULT  46"-66"

#### Lateral Supports (select one)

8"-16.5"W Range (select pad style below)

Flat Pads  Curved Pads

11"-19.5"W Range (select pad style below)

Flat Pads DEFAULT  Curved Pads

### PK207 Maximum Support Package.....\$7,989

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front for Shadow Tray, Adjustable/Removable Actuator Handle, Shadow Tray, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Positioning Belt, Accessories Mounting Bracket, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Tray, Multi-Adjustable Foot Plates, Planar Seat, Front Wheels/Casters, 5" Rear Locking Casters.

#### Front Wheel/Casters (select one)

5" Front Wheels DEFAULT  Front Swivel Casters

#### Actuator Handle Location (select one)

Right DEFAULT  Left

#### Shadow Tray (select one)

Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)

Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

#### Chest Pad (select one)

Non-Adjustable DEFAULT  Adjustable

#### Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT  21"L

#### Hip Supports (select one range)

10"-15"  12"-17" DEFAULT  17"-23"

#### Removable Back (select one)

Flat Back-15"H  Flat Back-19"H

Contoured Back-15"H  Contoured Back-19"H DEFAULT

#### Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT  38"-67"

#### Chest Support (select one)

Chest Strap select one chest measurement:  30"-46"  34"-54"  46"-66"

X-Style 11"Lx9.5"W  X-Style 14.5"Lx10.5"W DEFAULT

Y-Style 17.5"Lx11.5"W  Y-Style 20.5"Lx12.5"W

Non-Stretch X-Style 11"Lx9.5"W  Non-Stretch X-Style 14.5"Lx10.5"W

Non-Stretch Y-Style 17.5"Lx11.5"W  Non-Stretch Y-Style 20.5"Lx12.5"W

#### Lateral Supports (select one)

8"-16.5"W Range (select pad style below)

Flat Pads  Curved Pads

11"-19.5"W Range (select pad style below)

Flat Pads DEFAULT  Curved Pads

#### Head Support (select one)

Short  Tall DEFAULT

### PK208 Minimum Support Package (Mobile)....\$6,403

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

#### Actuator Handle Location (select one)

Right DEFAULT  Left

#### Chest Pad (select one)

Non-Adjustable DEFAULT  Adjustable

#### Removable Back (select one)

Flat Back-15"H  Flat Back-19"H DEFAULT

Contoured Back-15"H  Contoured Back-19"H

#### Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT  38"-67"

#### Chest Strap (select one, chest measurement)

30"-46"  34"-54" DEFAULT  46"-66"