

**Bill To** (Ultimate Medical Supplier Only)

|                                |
|--------------------------------|
| Purchase Order #:              |
| Quote #:                       |
| Account #:                     |
| Ordered By:                    |
| Supplier Company:              |
| Address:                       |
| City, State, Zip:              |
| Phone: _____ Email: _____      |
| Payer Source:                  |
| <b>Name of ATP (Required):</b> |

**Ship To**

|                   |
|-------------------|
| Name:             |
| Address:          |
| City, State, Zip: |
| Phone:            |
| Tag For:          |

**Glider Medium User Specifications**

|  |             |
|--|-------------|
| Approximate Height Range                           | 4'0"-5'6"   |
| Weight Limit                                       | 200 lbs.    |
| Seat Depth Range (center of knee pad to flat back) | 20"-25"     |
| Seat to Foot Plate Range (from seat pivot)         | 13.5"-16.5" |

**Base** *Note: The DEFAULT items will be added if no selection is made.*

|  |         |
|--|---------|
| <b>EasyStand<sup>®</sup> Glider Medium</b> | \$7,688 |
| <b>Price Includes:</b>                     |         |
| PNG50378 Glider Base                       | \$3,515 |
| PNG50187 Glider Configuration              | \$4,173 |

**Lift Mechanisms-must select one**

|  |         |
|--|---------|
| <input type="checkbox"/> PNG50113 Manual Hydraulic Actuator <b>DEFAULT</b> | N/C     |
| <b>Manual Hydraulic Actuator Handle Location</b>                           |         |
| <input type="checkbox"/> PNG50059 Right Side <b>DEFAULT</b>                | N/C     |
| <input type="checkbox"/> PNG50126 Left Side                                | N/C     |
| <b>Manual Hydraulic Actuator Handle Accessories</b>                        |         |
| <input type="checkbox"/> PNG50328 Adjustable/Removable Actuator Handle     | \$247   |
| <input type="checkbox"/> PNG50044 Quad Grip Handle Extension               | \$267   |
| <input type="checkbox"/> PNG50070 T-Style Handle Extension                 | \$167   |
| <input type="checkbox"/> PNG50010 Pow'r Up Lift (battery powered)          | \$1,364 |

**Tray-must select one**

|  |       |
|--|-------|
| <input type="checkbox"/> PNG50370 Glider Tray <b>DEFAULT</b> | \$362 |
|--|-------|

**Foot Plate Accessories**

|   |       |
|---|-------|
| <input type="checkbox"/> PNG30030 Foot Straps-15"L (length over top of foot, pair)            | \$84  |
| <input type="checkbox"/> PNG50505 Foot Straps-21"L (length over top of foot, pair)            | \$95  |
| <input type="checkbox"/> PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair) | \$136 |
| <input type="checkbox"/> PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair) | \$147 |

**Seats-must select one**

|   |       |
|---|-------|
| <input type="checkbox"/> PNG50350-1 Full Range Seat <b>DEFAULT</b>              | \$634 |
| <input type="checkbox"/> PNG50351 Full Range Seat with Roho <sup>®</sup> Insert | \$823 |

**Seat Accessories**

|   |       |
|---|-------|
| <input type="checkbox"/> PNG50499 Easy-Adjust Seat Depth- <b>CANNOT be ordered off-unit</b><br>(roller bearing allows for smoother seat depth adjustment) | \$118 |
|---|-------|

**Glide Handles**

|  |       |
|--|-------|
| <input type="checkbox"/> PNG50251 Glide Handle Extensions (replaces DEFAULT glide handles, pair) | \$183 |
|--|-------|

**Pelvic Support**

|  |       |
|--|-------|
| <input type="checkbox"/> PNG50345 Hip Supports (9"-15" range, pair)                                  | \$393 |
| <input type="checkbox"/> PNG50374 Hip Supports (11"-17" range, pair)                                 | \$393 |
| <input type="checkbox"/> PNG30029 Velcro <sup>®</sup> Positioning Belt-26"-55" (back required)       | \$84  |
| <input type="checkbox"/> PNG30306 Velcro <sup>®</sup> Positioning Belt-38"-67" (back required)       | \$95  |
| <input type="checkbox"/> PNG30028 Positioning Belt with Airline Style Buckle-26"-55" (back required) | \$183 |

**Backs**

|   |       |
|---|-------|
| <input type="checkbox"/> PNG50066 Flat Back-15"H                | \$721 |
| <input type="checkbox"/> PNG50067 Flat Back-19"H                | \$760 |
| <input type="checkbox"/> PNG50062 Contoured Back-15"H           | \$793 |
| <input type="checkbox"/> PNG50068 Contoured Back-19"H           | \$832 |
| <input type="checkbox"/> PNG50171 Removable Flat Back-15"H      | \$826 |
| <input type="checkbox"/> PNG50173 Removable Flat Back-19"H      | \$866 |
| <input type="checkbox"/> PNG50172 Removable Contoured Back-15"H | \$899 |
| <input type="checkbox"/> PNG50174 Removable Contoured Back-19"H | \$940 |

**Back Accessories**

|  |       |
|--|-------|
| <input type="checkbox"/> P600002 Back Angle Adjustment Wheel (back required) | \$162 |
| <input type="checkbox"/> P80246 Angle Locator                                | \$42  |

**Chest Straps-back required**

|   |       |
|---|-------|
| <input type="checkbox"/> PNG50192 Chest Strap-30"-46" (back required) | \$188 |
| <input type="checkbox"/> PNG50193 Chest Strap-34"-54" (back required) | \$198 |
| <input type="checkbox"/> PNG50504 Chest Strap-46"-66" (back required) | \$208 |

### Upper Body Support-back required

**ATTENTION - ALL OPTIONS BELOW REQUIRE THE ACCESSORIES MOUNTING BRACKET**

- PNG50034 Accessories Mounting Bracket \$100

### Chest Vests-back and accessories mounting bracket required

- PNG11212 High Mount Chest Vest Bracket \$85
- PNG30170 X-Style Chest Vest-11"Lx9.5"W \$321
- PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$327
- PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$321
- PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$327
- PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W \$371
- PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W \$377
- PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W \$371
- PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W \$377

### Lateral Supports-back and accessories mounting bracket required

- PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair) \$393
- PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair) \$393
- PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair) \$393
- PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair) \$393

### Head Supports-back and accessories mounting bracket required

- PNG50240 Head Support-Short \$427
- PNG50037 Head Support-Tall \$437

### Other Accessories

- PNG50040 Push Handles (back required, pair) \$167
- P80809 Tool Pouch \$29

### Notes:

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

**(A) Seat Depth Range: 20"-25"**  
(center of knee pad to flat back)

Measure the distance between the back of the buttocks and the back of the knee.

**(B) Seat to Foot Plate Range: 13.5"-16.5"**  
(from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



# AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

info@easystand.com easystand.com

Email form to [webquotes@easystand.com](mailto:webquotes@easystand.com)