



# Retail Order Form

USA Effective January 1, 2024 Date: \_\_\_\_\_

## Glider User Specifications

This price list is for both the Glider Medium and the Glider Large. Specifications for each are listed above the pricing.

### Bill To (EasyStand Supplier Only)

Purchase Order #: \_\_\_\_\_

Quote #: \_\_\_\_\_

Account #: \_\_\_\_\_

Ordered By: \_\_\_\_\_

Name of ATP (Required): \_\_\_\_\_

Supplier Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Payer Source: \_\_\_\_\_

### Ship To

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City, State: \_\_\_\_\_

Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Tag For: \_\_\_\_\_

Email form to [webquotes@easystand.com](mailto:webquotes@easystand.com)

# AltimateMedical

Ph: 800.342.8968 or 507.697.6393  
 Fax: 877.342.8968 or 507.697.6900  
[easystand.com](http://easystand.com)

# Glider Medium & Large Packages

## Glider Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (center of knee pad to flat back)	20"-25"
Seat to Foot Plate Range (from seat pivot)	13.5"-16.5"

## Glider Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (center of knee pad to flat back)	24"-29"
Seat to Foot Plate Range (from seat pivot)	15.5"-18.5"

## EasyStand Glider Medium

**PK213 Moderate Support Package.....\$10,859**

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Glide Handle Extensions, Removable Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Adjustable Foot Plates.

### Actuator Handle Location (select one)

- Right DEFAULT
- Left

### Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT
- 21"L

### Hip Supports (select one range)

- 9"-15" DEFAULT
- 11"-17"

### Removable Back (select one)

- Flat Back-15"H
- Flat Back-19"H DEFAULT
- Contoured Back-15"H
- Contoured Back-19"H

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT
- 38"-67"

### Chest Strap (select one, chest measurement)

- 30"-46"
- 34"-54" DEFAULT
- 46"-66"

## EasyStand Glider Large

**PK214 Moderate Support Package.....\$10,687**

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Removable Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Glide Handles, Adjustable Foot Plates.

### Actuator Handle Location (select one)

- Right DEFAULT
- Left

### Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT
- 21"L

### Hip Supports (select one range)

- 10"-15"
- 12"-17" DEFAULT

### Removable Back (select one)

- Flat Back-15"H
- Flat Back-19"H DEFAULT
- Contoured Back-15"H
- Contoured Back-19"H

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55"
- 38"-67" DEFAULT

### Chest Strap (select one, chest measurement)

- 30"-46"
- 34"-54" DEFAULT
- 46"-66"