

## MSRP Itemized Price List/Order Form

USA Effective July 1, 2019 Date: \_\_\_\_\_

### Bantam Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat pivot to planar back)	15"-20"
Seat to Foot Plate Range (from seat pivot)	7"-18.5"

### Bill To (Altimate Medical Supplier Only)

Purchase Order #: \_\_\_\_\_

Quote #: \_\_\_\_\_

Account #: \_\_\_\_\_

Ordered By: \_\_\_\_\_

Name of ATP (Required): \_\_\_\_\_

Supplier Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Payer Source: \_\_\_\_\_

### Ship To

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City, State: \_\_\_\_\_

Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Tag For: \_\_\_\_\_

**Email form to [webquotes@easystand.com](mailto:webquotes@easystand.com)**

# Itemized Price List

PY5500-1 EasyStand<sup>®</sup> Bantam Medium.....\$1,998

### The default includes:

- Two 5" Front Wheels
- Manual Hydraulic Actuator with Handle
- Height-Only Adjustable Foot Plates
- Blue Frame Color
- Tray and Back Support Required
- Two 5" Rear Locking Casters
- Planar Seat
- Black Upholstery

### Choose actuator handle location:

- Right Side (default)  
 Left Side

### Foot Holder Size (please select one)

- 7.75"Lx3.25"W  
 9.75"Lx4"W  
 11.75"Lx5"W (default)

### Knee Pad Size (inside width, please select one)

- 4.0"  
 5.5" (default)

### Supine (CANNOT be ordered off-unit)\*

- PY5554 Supine \$514  
(shadow tray, back option and head support required)

### Seat Adjustment (CANNOT be ordered off-unit)

- PY5645 Easy-Adjust Seat Depth \$98  
(Roller bearing seat depth that allows for smoother seat depth adjustment)

### Front Casters

- PY5568 Front Swivel Casters \$163

### Lift Mechanisms

- PY5564 Pow'r Up Lift (replaces hydraulic actuator, battery powered) \$1,119  
(not available with pow'r up)  
 PY5632 Adjustable/Removable Actuator Handle \$208  
 PNG50044 Quad Grip Handle Extension (not available with pow'r up) \$140  
 PNG50070 T-Style Handle Extension (not available with pow'r up) \$140

### Shadow Trays (must select one)

- PY5604 Black Molded Tray-21"Lx24"W (depth, height and angle adjustable) \$630  
 PY5606 Clear Tray-21"Lx24"W (depth, height and angle adjustable) \$871  
 PY5656 Black Molded Swing-Away Tray-21"Lx24"W (depth, height and angle adjustable) \$747  
 PY5657 Clear Swing-Away Tray-21"Lx24"W (depth, height and angle adjustable) \$922

### Tray Accessories

- PNG50042 Hand Grips (for clear tray only, pair) \$140  
 PT50156 Elbow Pad-9.25" Cutout \$131  
 PT50160 Elbow Pad-11" Cutout \$131  
 PY5582 Elbow Pad-13" Cutout \$131  
 PNG30000 Large Contoured Chest Pad (not available with elbow pads) \$129  
 PNG30295 Padded Tray Cover (not available with elbow pads) \$206

### Knee Pads-must select one

- PY5594-1 Slide-In Knee Pads (if no selection is made, slide-in knee will be selected, pair) \$400  
 PY5596-1 Swing-Away Knee Pads (pair) \$644

### Knee Accessories

- PY5600-1 Roho<sup>®</sup> Knee Pads (pair, only available with 5.5" pad) \$285

### Foot Plates and Foot Plate Accessories

- PY5586 Multi-Adjustable Foot Plates (pair) \$229  
 PY5636 Foot Straps-15"L (length over top of foot, pair) \$41  
 PY5638 Secure Foot Straps-15"L (length over top of foot, two pair) \$80

### Seats and Seat Accessories

- PY5570 Contoured Seat-Narrow 11"-13"W \$369  
 PY5572 Contoured Seat-Wide 12"-14"W \$369  
 PY3042 Black Hygienic Seat Cover for Planar Seat \$85  
 PY3043 Black Hygienic Seat Cover for PY5570 \$85  
 PY3044 Black Hygienic Seat Cover for PY5572 \$85  
 PY5560 Hip Supports (9"-15" range, pair, not available with contoured seat) \$332

### Backs (must select a one)

- PY5530 Planar Back-16"H \$696  
 PY5532 Planar Back-20"H \$735  
 PY5574 Contoured Back-16"H \$767  
 PY5576 Contoured Back-20"H \$806

### Back Accessories

- PY3041 Black Hygienic Back Cover for PY5530 \$85  
 PY3040 Black Hygienic Back Cover for PY5532 \$85  
 PY3045 Black Hygienic Back Cover for PY5574 \$85  
 PY3046 Black Hygienic Back Cover for PY5576 \$85

### Positioning and Upper Body Support

- PY3008 Velcro<sup>®</sup> Positioning Belt (fits hip circumference up to 66", velcro with D-ring) \$51  
 PY3014 Positioning Belt with Airline Style Buckle (fits hip circumference up to 66", velcro with D-ring) \$140  
 PY5620 Chest Strap (circumference up to 42", velcro with D-ring) \$129  
 PY5634 High Mount Chest Vest Bracket (head support required) \$71  
 PNG30170 X-Style Chest Vest-11"Lx9.5"W \$253  
 PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$253  
 PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$253  
 PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$253  
 PY5562 Lateral Supports-9"-14.5" Range (pair, not available with contoured back) \$332  
 PY5658 Lateral Supports with Elbow Stop-9"-14.5" Range (for planar pads only, pair) \$339  
 PY5659 Lateral Supports with Elbow Stop and Arm Rest-9"-14.5" Range (for planar pads only, pair) \$375  
 PY5660 Elbow Stop with Arm Rest (for planar and contoured pads, pair) \$339

### Head Supports

- PY5626 Head Support-5"Hx8"W \$332  
 PY5628 Head Support-6"Hx10"W \$332  
 PY5630 Head Support-7"Hx14"W \$332  
 PA5586 Head Support-Form to Fit-5"Hx10"W \$332

### Other Accessories

- PNG50040 Push Handles (pair) \$140  
 P80809 Tool Pouch \$22  
 PY5643 Angle Locator \$32

# AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

[easystand.com](http://easystand.com)